SMARTBMI SCALE BT

MT5513


User manual

EN

Introduction


Thank you for buying this modern scale with very interesting features. Device can work as simple scale, or smart scale when you use with smartphone – device will connect via Bluetooth and send data directly to smartphone (please download and install free application “OKOK” first). Application can show you BMI, body fat factor and many more. Please read this manual before using the device.

Overview

Getting started

1. Install “OKOK” application (it is a free application available in GooglePlay™ for Android smartphone).

2. Enable Bluetooth interface in your smartphone and stand on the scale with bare feet.
3. Scale will automatically connect with application.
4. Wait until for measurement is done.

Note! Scale lets a small amount of weak (less than 1mA) current flow through the human body, so as to detect the bioimpedance and estimate body fat. However, please be aware that anyone with an wearable or implantable medical electronic instrument, **such as a pacemaker**, must avoid using this device.

Due to continuous development specification and appearance of product are subject to change without prior notice. For technical support please visit www.media-tech.eu.

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EN Environmental guidance:

The device is marked with the symbol of a crossed out waste container, in accordance with the European Directive 2012/19/EC on waste electrical and electronic equipment. Products marked with this symbol should not be disposed of or disposed of after the end of their useful life together with other household waste. The user is obliged to dispose of waste electrical and electronic equipment by delivering it to a designated point where such hazardous waste is recycled. The collection of this type of waste in separate places and the proper process of their recovery contribute to the protection of natural resources. Proper recycling of waste electrical and electronic equipment has a beneficial effect on human health and the environment. In order to obtain information on where and how to dispose of waste electrical and electronic equipment in an environmentally safe way, the user should contact the relevant local authority, the waste collection point or the point of sale where he purchased the equipment.

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